## Flow Focusing Monodisperse Aerosol Generator

## **Model 1520**

Produces monodisperse droplets and solid particles.

The Flow Focusing Monodisperse Aerosol Generator (FMAG) 1520 uses the aerodynamic flow-focusing effect to accurately control the diameter of a liquid jet for generating monodisperse droplets from 15 to 90 µm in diameter, which can then be dried to produce particles from 0.8 to 12 µm in diameter. In normal operation, a built-in syringe pump pushes liquid out of a 100-µm-diameter nozzle and is stretched to a much thinner stream by the focusing gas flow. The resulting thin liquid jet then breaks up into uniform-sized droplets after passing through a vibrating ceramic aerosol generation head. A coaxial flow of clean air is introduced to dry the droplets into solid particles. The aerosol then exits the top of the FMAG after passing by a built-in electrical corona ionizer and inspection light for easy viewing. The large 100-µm-diameter nozzle in the FMAG enables aerosol generation over extended periods of time without experiencing nozzle clogging problems, and at a very low liquid pressure. This low shear stress generally enables biological cells to remain viable, even after dispersion as uniform particles. The closely-related model 1530 omits the drying air, ionizer, and light, and has provisions for remotely mounting the aerosol generation head upside down for generating large droplets.



## **Polydisperse Generators**

This type of generator is typically capable of spraying aqueous solutions (e.g. salt), suspensions PSL, gold or glass nanoparticles), or oil or similar substances. They serve a variety of applications from laboratory research, field tests of detectors, and filter testing.

Model	3073	3079A	3076	8026	8108
Particle Size Range (µm)	0.01 to 2.0 (nominal 0.3 count mean diameter)				0.1 to 10
Particle Concentration (particles/cm³) in Output	10 <sup>2</sup> to 10 <sup>7</sup>	>108	>107	10 <sup>6</sup>	<10³ at 1 µm
Nominal Flow Rate (L/min)	0.3 to 4.5	1.0 to 4.2	3.0	1.25	140
Note	Portable, battery option available	Portab <b>l</b> e	Laboratory grade	Portab <b>l</b> e	Designed for ISO 16890-2 and ASHRAE 52.2 filter testing

